


# Rumba & Rumba modular








1BA	RUMBA MODULAR
1071	NOME/ (MODUL) III

SEATING	4	SEATING	14
BACKREST /		BACKREST	14
BACKREST RECLINE OPTIONS	4		
		BACKREST RECLINE OPTIONS	15
ARMREST	5		
		ARMREST	15
LEGRESTS			
Swing Away & Elevating Legrests	5	JAY CUSHIONS & BACKRESTS	16
Footplates	5		
		LEGRESTS	
JAY CUSHIONS & BACKRESTS	6	Swing Away & Elevating Legrests	18
		Footplates	18
CONTROL SYSTEMS	8		
VR2 Controls	8	CONTROL SYSTEMS	19
Control Mountings	9	VR2 Controls	19
		PG R-Net Controls	21
WHEELS	10	Control Mountings	23
BATTERIES & CHARGERS	11	WHEELS	24
ACCESSORIES	11	BATTERIES & CHARGERS	25
		ACCESSORIES	26
		Lights & Indicators	27

# Seating







# Backrest/BackrestReclineOptions



Standard sling backrest



Tension adj. sling backrest



Back angle ajustable with Tools

The standard rehab manual adjustable back provides a fixed back angle from 0° to 30° Adjustment is still quick and easy. Simply remove the fixing bolts, set the desired angle and refit the bolts.



## **Armrest**



Std Armrest (height Adjustable) (+ 3cm width adj. each side)

# Legrests

## Swing Away & Elevating Legrests



Swing-away hangers

Push the button to swing away and detach. The range of swing away legrests is available in all frame colours.

70° hanger



Manual elevating legrest

The manual elevating legrest can be easily adjusted by the flip of a lever. Swing-away, detachable and complete with knee pads as standard.

## **Footplates**



Composite footplates standard for Rumba

# Jay clinical cushions and backs

JAY° offers a complete range of cushions and backs to meet individual user needs. The comfort cushions, for users with moderate to low skin protection and positioning needs. Skin protection cushions for moderate sking protection and positioning needs and the Positioning & Skin ptotection cushions, for user with more complex positioning needs and or high skin breakdown risk

For more information please refer to the JAY® cushions Brochure and Order form





#### **Skin Protection**



### Positioning & Skin protection, Positioning






# Jay Backs

JAY° backs are designed to provide posterior lateral pelvic stabilization along with enhanced trunk stability. They are lightweight and adjustable, allowing for easy installation. JAY° backs reduce fatigue and maximize function to provide clients the comfort and

stability to go anywhere with confidence.

For more information please refer to the JAY\* backs Brochure and Order form



# Control-Systems

## **VR2 Controls**

VR2 Remote control Attendant only (it is the same as Standard Remote control)



VR 2 with actuator controls





VR2 lights & indicators



VR2 Dual (User + Attendant)



VR2 Directional control (User + Attendant)

VR2 Directional control (attendant only))





## **Control Mountings**

## Fixed front mount

The fixed front mount is the standard control mount for Quickie hand controls.



## Swing-away front mount

Unlike the fixed front mount this swing-away mount offers clients horizontal adjustment of the hand control allowing them to swing it away and adjust it to their preference.



## Chin control



# Wheels

8,5" x 2" Solid 8,5" x 2" Pneumatic





9" x 2.5" Solid 9" x 2.5" Pneumatic





12.5" Solid

# **Batteries & Chargers**

50 Ah maintenance free 40 Ah maintenance free



220V / 8A automatic charger





# Accessories

Crutch holder



5 cm (2") positioning strap



Swing Away tray



Kerb climber with 8" castor wheels

Kerb climber with 9" castor wheels



# Accessories

## **Headrests**

Standard Headrest with sky upholstery with Anatomic backrest



Anatomic Headrest with sky upholstery with Anatomic backrest



# Lights & Indicators



# Rumba Modular



# Seating

5 cm (2") standard foam cushion Our standard foam cushion is available in sizes to suit all seat widths





## **Backrest**



Standard sling backrest



Tension adj. sling backrest (black)



## **Backrest Recline Options**



Back angle ajustable with Tools

The standard rehab manual adjustable back provides a fixed back angle from 0° to 30° Adjustment is still quick and easy. Simply remove the fixing bolts, set the desired angle and refit the bolts.



Manual recline – gas spring system Only available for Rumba Modula

If your client requires a back that can be reclined whilst in the chair, then the gas spring system is for them.

With the client seated in their Rumba, carers can recline this back from 0° to 30° with the press of a lever.

With stabilization bar, with Push handlers



Power recline Only available for Rumba Modula

Ideally suited for clients who want to adjust the back angle with no assistance from a carer.

Powered from the hand control, the back can be reclined from 0° to 30°.

With stabilization bar, with Push handlers

## **Armrest**

Std Armrest (height Adjustable)

(+ 3cm width adj. each side)





# Jay clinical cushions and backs

JAY® offers a complete range of cushions and backs to meet individual user needs. The comfort cushions, for users with moderate to low skin protection and positioning needs. Skin protection cushions for moderates king protection and positioning needs and the Positioning & Skin ptotection cushions, for user with more complex positioning needs and or high skin break down risk

 $For more information please refer to the JAY ^\circ cushions Brochure and Order for matter than the properties of the prop$ 





#### **Skin Protection**



### Positioning & Skin protection, Positioning






# Jay Backs

JAY° backs are designed to provide posterior lateral pelvic stabilization along with enhanced trunk stability. They are lightweight and adjustable, allowing for easy installation. JAY° backs reduce fatigue and maximize function to provide clients the

comfort and stability to go anywhere with confidence.

For more information please refer to the JAY° backs Brochure and Order form



# Legrests

## Swing Away & Elevating Legrests



Swing-away hangers

Push the button to swing away and detach. The range of swing away legrests is available in all frame colours.

70° hanger



Manual elevating legrest

The manual elevating legrest can be easily adjusted by the flip of a lever.

Swing-away, detachable and complete with knee pads as standard.



Powered elevating legrest (only Available in Rumba Modular)

The powered elevating legrest is powered from the hand control and gives even the most challenged client the ability to adjust the legrest.

Swing-away, detachable and complete with knee pads as standard.

## **Footplates**





Composite footplates Aluminium footplate



# ControlSystems

## **VR2 Controls**

VR2 Remote control Attendant only (it is the same as Standard Remote control)



VR 2 with actuator controls



VR2 with actuator controls and lights & indicators



VR2 lights & indicators



# Control-Systems

VR2 Dual (User + Attendant)



VR2 Directional control (User + Attendant)

VR2 Directional control (attendant only))







## **PG R-Net Controls**

R-Net LED Joystick, with Actuator Controls

Multiple profile hand control used for drive and actuator control. Assignable Button function built into this joystick.





R-Net LED Joystick, with Actuator Controls & Lights

Multiple profile hand control used for drive, actuators and lighting control. Assignable Button function built into this joystick.





#### **R-Net Dual Control**

Compact multiple profile hand control mounted at the rear for an attendent.

Control profiles and seating functions with this joystick.

R-Net only attendant



R-NET -Directional control (attendant only)

R-Net - Standard Module



R-Net Lighting Module (ISM-L)

This module is required when Lights & Indicators are desired.

R-Net Seating & Lighting Module (ISM-6L)

This seating and lighting module is required to upgrade the system from one to a maximum of six powered seating functions, whilst also providing lights and indicators function.





## **Control Mountings**

## Fixed front mount

The fixed front mount is the standard control mount for Quickie hand controls.



## Swing-away front mount

Unlike the fixed front mount this swing-away mount offers clients horizontal adjustment of the hand control allowing them to swing it away and adjust it to their preference.



# Wheels

## **Caster Wheels**



8,5" x 2" Solid 8,5" x 2" Pneur



9" x 2.5" Solid 9" x 2.5" Pneumatic

## **Drive Wheels**





12.5" Solid
12.5" Pneumatic
This slick tread is ideal for indoor and outdoor use. Available in solid and pneumatic.

# **Batteries & Chargers**

50 Ah maintenance free 40 Ah maintenance free



220V / 8A automatic charger





# Accessories

Crutch holder



5 cm (2") positioning strap



Swing Away tray



Kerb climber



# Accessories

## **Headrests**



Standard Headrest with sky upholstery with Anatomic backrest



Anatomic Headrest with sky upholstery with Anatomic backrest





# Lights & Indicators



















All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com or alternatively is available on request in large text.

Sunrise Medical GmbH & Co. KG

Sunrise Medical S.L Polígono Bakiola, 41 48498 Arrankudiaga Vizcaya - Spain

Tel: +34 (0) 902142434 www.quickiepower.com For international enquiries e-mail: International@sunrisemedical.de